


# Tubbing Guide

---



2 or more 15+ Qt/3.75g food-grade plastic tubs	\$7/each Target
Seachem Prime	\$10/8.5 floz Petco
1ml Graduated Transfer Pipettes	\$20/box of 500
Thermometer	\$5
Egg Crate or Plastic Canvas	Lowe's or Joann's
Appropriate Hide	-
Air stone & Air pump	-
Indian Almond Leaf	-

---

## Preparation:

I suggest using Sterilite 15 qt tubs from Target.  
Fill one of the tubs  $\frac{2}{3}$  to  $\frac{3}{4}$  full of 60-65°F water.  
Add 0.25mL of Seachem Prime using graduation on transfer pipette.  
Note that RO, DI, and soft water need additional chemical adjustment.  
Place tub in dark, cool location.  
Add hide, air stone, and air pump.  
Add IAL.

---

## Procedure:

Place axolotl in the tub. One axolotl per tub.

Adjust the flow of the air stone and air pump for axolotl's comfort.

Change 100% of the water every 12-24 hours. This is where having two tubs is helpful. I recommend preparing the new tub, then moving the axolotl over, then cleaning out the old tub.



Lost on a Mountain  
Axolotls




loamaxolotls@gmail.com



LostonaMountain.com

Check out my basic video at TikTok: <https://www.tiktok.com/@lostmountainaxolotl/video/7066958304140086574>



---

Tubbing can be used while you cycle your tank (see cycling guide), during quarantine (all new axolotls should be quarantined for 30 days before adding to a tank with another axolotl – same gender within 1” in size), or if you need to do health treatments.

